

# Serum Amyloid A (SAA) Screening for General Health and Biosecurity

Author: Siddra Hines, DVM, PhD, DACVIM-LA

**AUG**  
2020

## Relevant Species



Routine equine health exams and biosecurity strategies utilize screening methods to detect underlying performance-limiting issues and potential infectious disease concerns. The fundamental basis of this is a good physical exam with rectal temperature screening, however subclinical issues may still exist that can become problematic. Measurement of serum amyloid A (SAA) is a highly effective tool for this purpose and can identify horses with subclinical disease that would not otherwise be detected.



## Any elevation of SAA should be considered abnormal

The normal range of SAA in a healthy horse is generally considered to be less than 20  $\mu\text{g}/\text{mL}$ .<sup>1,2</sup> It becomes elevated in acute systemic inflammation, particularly due to viral or bacterial infections, with bacteria stimulating the greatest SAA production.<sup>3-5</sup> SAA is highly sensitive and specific for the presence or absence of systemic inflammation, more sensitive than WBC count or fibrinogen.<sup>1,3,6-8</sup> For the purpose of biosecurity or health screening, any elevation of SAA should be considered abnormal and trigger further investigation.<sup>2,7</sup>

## SAA may increase even if body temperature is normal

SAA can increase even in the absence of fever<sup>9-10</sup>, making it exceptionally useful for screening purposes and to monitor health status. Unlike fever, it will not be significantly affected by NSAID therapy.<sup>7</sup> In one study, SAA was 97.1% sensitive and 97.2% specific to differentiate clinically abnormal horses (i.e. those who developed infections) from those that were normal when tested 24 hours after air transportation.<sup>9</sup> The presence of fever had a sensitivity of just 2.9% at the same time point, indicating that SAA could identify brewing infections earlier and allow more rapid clinical intervention.<sup>9</sup>

## SAA testing can help prevent and manage infectious disease outbreaks

There is clear value in assessing SAA alongside rectal temperature to identify subclinical disease<sup>2,3,9</sup> as part of an efficacious health screening strategy. In terms of biosecurity, disease outbreaks may be prevented if horses are effectively screened prior to co-mingling at events or when introducing new horses to a resident population. If an active infectious disease outbreak is occurring, SAA can help monitor at-risk or exposed horses for development of disease. Horses that have become infected will have elevated SAA and can be handled appropriately, with additional diagnostics as indicated.<sup>3,4</sup> SAA can also be used to monitor populations that may be at increased risk due to age, stress, exposure, population density, or other factors.<sup>4,8</sup> This could include young horses in intense training<sup>4</sup>, hospital populations<sup>6</sup>, or horses undergoing long-distance travel<sup>9</sup>.

## SAA can identify subclinical problems in outwardly healthy horses

As a general health screening tool, SAA can be very useful to assess the health of a horse prior to surgery and to monitor for complications afterwards.<sup>3,10-12</sup> Testing prior to transport may identify subtle abnormalities that have the potential to develop into bigger problems with shipping. Testing horses after transport<sup>3,9</sup> allows early recognition of shipping-related infections. SAA screening can help detect underlying issues during any routine health examination, including pre-purchase or insurance exams. Although normal (negative) results do not rule out all concerns, a positive result indicates an active problem that should be investigated.

## Normal horses should have little to no SAA

As a general rule, a truly normal horse should have virtually no circulating SAA. SAA increases minimally or not at all with stress, exercise<sup>13,14</sup> or anesthesia<sup>4</sup> alone. Barring any confounding factors (see April 2020 newsletter), elevated SAA in an outwardly healthy horse should always be a trigger to look deeper into possible causes.

---

## REFERENCES

<sup>1</sup> Witkowska- Piłaszewicz OD, Żmigrodzka M, Winnicka A, et al. Serum amyloid A in equine health and disease. *Equine Vet J* 2019;51(3):293-298.

<sup>2</sup> Radcliffe RM, Buchanan BR, Cook VL, et al. The clinical value of whole blood point-of-care biomarkers in large animal emergency and critical care medicine. *J Vet Emerg Crit Care* 2015;25(1):138-151.

<sup>3</sup> Jacobsen S, Andersen PH. The acute phase protein serum amyloid A (SAA) as a marker of inflammation in horses. *Equine Vet Educ* 2007;19(1):38-46.

<sup>4</sup> Pepys MB, Baltz ML, Tennet GA. Serum amyloid A protein (SAA) in horses: objective measurement of the acute phase response. *Equine Vet J* 1989;21(2):106-109.

<sup>5</sup> Hultén C, Demmers S. Serum amyloid A (SAA) as an aid in the management of infectious disease in the foal: comparison with total leucocyte count, neutrophil count and fibrinogen. *Equine Vet J* 2002;34(7):693-698.

<sup>6</sup> Belgrave RL, Dickey MM, Arheart KL. Assessment of serum amyloid A testing of horses and its clinical application in a specialized equine practice. *J Am Vet Med Assoc* 2013;243(1):113-119

<sup>7</sup> Hooijberg EH, van der Hoven R, Tichy A, et al. Diagnostic and predictive capability of routine laboratory tests for the diagnosis and staging of equine inflammatory disease. *J Vet Intern Med* 2014;38:1587-1593.

<sup>8</sup> Anhold H, Candon R, Chan Di-Sien, et al. A comparison of elevated blood parameter values in a population of thoroughbred racehorses. *J Equine Vet Sci* 2014;34:651-655.

<sup>9</sup> Oertly M, Gerber V, Anhold H, et al. The accuracy of serum amyloid A in determining early inflammation in horses following long-distance transportation by air, in Proceedings. *Am Assoc Equine Pract* 2017;63:460-461.

<sup>10</sup> Jacobsen S, Jensen JC, Frei S, et al. Use of serum amyloid A and other acute phase reactants to monitor the inflammatory response after castration in horses: a field study. *Equine Vet J* 2005;37(6):552-556.

<sup>11</sup> De Cozar M, Sherlock C, Knowles E, et al. Serum amyloid A and plasma fibrinogen concentrations in horses following emergency exploratory celiotomy. *Equine Vet J* 2020;52(1):59-66.

<sup>12</sup> Aitken MR, Stefanovski D, Southwood LL. Serum amyloid A concentration in postoperative colic horses and its association with postoperative complications. *Vet Surg* 2018; 1-9. <https://doi.org/10.1111/vsu.13133>

<sup>13</sup> Kristensen L, Buhl R, Nostell K, et al. Acute exercise does not induce an acute phase response (APR) in Standardbred trotters. *Can J Vet Res* 2014;78(2):97-102.

<sup>14</sup> Cywinska A, Witkowski L, Szarska E, et al. Serum amyloid A (SAA) concentration after training sessions in Arabian race and endurance horses. *BMC Vet Res* 2013;9:91-97.