

Elevate your patient care with vmrdSAA



In-Clinic Test

for rapid, objective assessment



Numerical results

with handheld reader


Serum amyloid A (SAA) is an acute phase protein that increases with acute systemic inflammation. As a biomarker, it is particularly valuable to help diagnose infection and other severe inflammatory diseases, by narrowing your differential list and objectively monitoring efficacy of subsequent treatment. Currently used in horses and cats, veterinarians can now utilize this valuable tool in their canine patients with the VMRD SAA product for Canine!

Multiple tests, one **handheld** reader - now with **RECHARGEABLE** battery!

- ✓ Ideal size with convenient storage case
- ✓ Calibration for accuracy & consistency
- ✓ Will not lose your results
- ✓ Automatically reads test results at 10 minutes



Cats can be hard to read. **vmrdSAA** can help.



In-Clinic Test

for rapid, objective assessment



Numerical results

with handheld reader

Serum amyloid A (SAA) is an acute inflammatory biomarker that can help diagnose infection and monitor treatment efficacy. Used extensively in horses for years, veterinarians can now utilize this valuable tool in cats with the VMRD SAA product for Feline!

Multiple tests, one handheld reader - now with RECHARGEABLE battery!

- ✓ Ideal size with convenient storage case
- ✓ Calibration for accuracy & consistency
- ✓ Will not lose your results
- ✓ Automatically reads test results at 10 minutes

