

**CERTIFICATE OF ANALYSIS**

# A3B10


Monoclonal Antibody

Catalog No.:	A3B10
Lot:	P231207
Expiration:	1 year when stored at 2-7 °C

**Quality Control Method:** IFA using CPV 12-well slide (catalog no. SLD-IFA-CPV), Isotype control IgG<sub>2a</sub> (1.0 mg/ml) and Anti-Mouse affinity purified FITC conjugate (catalog no. CJ-F-MURG-AP-10ML). RID using in-house IgG<sub>2a</sub> plate.

**Specific Reaction:** 3-4+ signal at 10 µg/ml with an endpoint of 0.1 µg/ml. The isotype control was negative with no background. The RID test showed a concentration of 1.097 mg/ml.

**Other Comments:** NA

**QC Manager:** 

David Hyatt

**Pattern Of Fluorescence:** Single cells with cytoplasmic and nuclear fluorescence.

Technical Data Sheet Version: 1