

CERTIFICATE OF ANALYSIS

Feline Leukemia Virus (FeLV)

Primary Antibody

Catalog No.:	AB1-FELV
Lot:	P170103-001
Expiration:	06 January 2021

Quality Control Method: Indirect FA using FeLV 2-well slide (catalog no. SLD-FAC-FELV), positive

> blood smear (catalog no. SLD-BSP-FELV), three negative blood smears (catalog no. SLD-BSN-FELV), and FeLV secondary reagent (catalog no.

AB2-FELV).

Specific Reaction: 1-2+ signal on positive smear, negative on the negative smear with no

background. 1-3+ signal on the positive well and negative on the negative

well with no background.

Other Comments: NA

Will Thurman QC Manager:

QC Release Date: 09 January 2017

Pattern Of Fluorescence: On the 2-well slide there should be Individual cell cytoplasmic fluorescence

with occasional small syncytia. On positive blood smear there is signal on

white blood cells and platelets (some clumped).

FA Cross-Reaction This product has been tested by FA and was found to react with Bartonella

Screening: Henselae at a 1+ signal and with Toxoplasma gondii at a 1+ signal but

does not react with Feline Calicivirus (FCV), Feline Herpesvirus (FHV-1/FVR), Feline Immunodeficiency Virus (FIV), Feline Infectious Peritonitis type 1 and 2 (FIP-1 and 2), Feline Panleukopenia Virus (FPLV), Reovirus (REO), and Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis

Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Conjugate\FeLV\FeLV Primary Antiserum\AB1-FELV_P170103-001_210106.docx

11 January 2017