

**CERTIFICATE OF ANALYSIS**

## Feline Leukemia Virus (FeLV)


Secondary FITC Conjugate

Catalog No.:	AB2-FELV
Lot:	P150409-001
Expiration:	21 April 2019

**Quality Control Method:** Indirect FA using FeLV 2-well control slide (catalog no. SLD-FAC-FELV), positive blood smear (catalog no. SLD-BSP-FELV), negative blood smear (catalog no. SLD-BSN-FELV), and primary anti-FeLV antibody (catalog no. AB1-FELV).

**Specific Reaction:** 2-4+ fluorescence with trace background on the positive well and negative with trace background on the negative well. 2-3+ fluorescence, with no background on the positive blood smear.

**Other Comments:** NA

**QC Manager:**  Torrey Beeson

**QC Release Date:** 22 April 2015

**Pattern Of Fluorescence:** On the 2 well slide there should be Individual cell cytoplasmic fluorescence with occasional small syncytia.

**FA Cross-Reaction Screening:** The raw material has also been tested satisfactorily by VMRD's IFA systems for *Bartonella henselae*, Feline Calicivirus (FCV), Feline Leukemia Virus (FeLV), Feline Herpesvirus (FHV-1/FVR), Feline Immunodeficiency Virus (FIV), Feline Infectious Peritonitis type 1 and 2 (FIP-1 and 2), Feline Panleukopenia Virus (FPLV), Reovirus (REO), and Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 1

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Quality VMRD\COAs and TDSCOA (VMRD)\Conjugate\FeLV\FELV Secondary Conjugate\AB2-FELV\_P150409-001\_190421.docx

23 April 2015

**FOR *IN VITRO* LABORATORY USE ONLY.**

H:\Quality VMRD\COAs and TDSCOA (VMRD)\Conjugate\FeLV\FeLV Secondary Conjugate\AB2-FELV\_P150409-001\_190421.docx

23 April 2015