

CERTIFICATE OF ANALYSIS

## Feline Leukemia Virus (FeLV)

Secondary FITC Conjugate

Catalog No.:	AB2-FELV
Lot:	P161215-004
Expiration:	15 July 2020

**Quality Control Method:** Indirect FA using FeLV 2-well control slide (catalog no. SLD-FAC-FELV), positive blood smear (catalog no. SLD-BSP-FELV), negative blood smear (catalog no. SLD-BSN-FELV), and primary anti-FeLV antibody (catalog no. AB1-FELV).

**Specific Reaction:** 2-4+ signal on the positive well, negative on the negative well with trace background on both wells. 1-2+ signal on the positive blood smear, negative on the negative blood smear.

**Other Comments:** NA

**QC Manager:**  Will Thurman

**QC Release Date:** 25 July 2016

**Pattern of Fluorescence:** On the 2 well slide there should be Individual cell cytoplasmic fluorescence with occasional small syncytia.

**FA Cross-Reaction Screening:** The raw material has also been tested satisfactorily by VMRD's IFA systems for Bartonella henselae, Feline Calicivirus (FCV), Feline Leukemia Virus (FeLV), Feline Herpesvirus (FHV-1/FVR), Feline Immunodeficiency Virus (FIV), Feline Infectious Peritonitis type 1 and 2 (FIP-1 and 2), Feline Panleukopenia Virus (FPLV), Reovirus (REO), and Rabies Recombinant Nucleoprotein (rNP), Toxoplasma gondii, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Conjugate\FeLV\FeLV Secondary Conjugate\AB2-FELV\_P161215-004-200715.docx

19 December 2016