

CERTIFICATE OF ANALYSIS

BA-26(A)

Monoclonal Antibody

Catalog No.:	BA-26(a)
Lot:	P240226
Expiration:	1 year when stored at 2-7 °C

Quality Control Method: Indirect FA using BVDV-1(Modderman), BVDV-2 (TN131) 12-well (in-house slides), isotype control IgG₁ and anti-murine affinity purified FITC conjugate (catalog no. CJ-F-MURG-AP-1ML or 10ML). RID using in-house IgG₁ plate.

Specific Reaction: Both BVDV type 1 and type 2 produced 2–3+ fluorescence at 10 µg/ml with an endpoint concentration below 1 µg/ml. The isotype control was negative with trace background. The RID test results showed a concentration of 1.024 mg/ml.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

21 March 2024