

CERTIFICATE OF ANALYSIS

Bovine Adenovirus Type 3 (BAV-3)

FITC Conjugate

Catalog No.:	CJ-F-BAV3-10ML
Lot:	P180502-004
Expiration:	10 May 2022

Quality Control Method: Direct FA using BAV-3 2-well slide (catalog no. SLD-FAC-BAV3).

Specific Reaction: 2 – 4+ signal on the positive well and negative on the negative well, trace background.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 10 May 2018

Pattern Of Fluorescence: Individual cells with diffuse cytoplasmic and membrane fluorescence and strong nuclear fluorescence.

FA Cross-Reaction Screening: The ready to use conjugate has also been screened by FA and has been found to react with BAV-1 at 1-2+ fluorescence, but does not react with Bovine Adenovirus Type 5 (BAV-5), Bovine Coronavirus (BCV), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Bovine Herpesvirus type 1 (BHV-1/IBR), Parainfluenza Virus Type 3 (PI-3), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

23 May 2018