

CERTIFICATE OF ANALYSIS

Bovine Adenovirus Type 3 (BAV-3)

FITC Conjugate

Catalog No.:	CJ-F-BAV3-10ML
Lot:	P230402-001
Expiration:	28 March 2031

Quality Control Method: Direct FA using BAV-3 2-well slide (catalog no. SLD-FAC-BAV3).

Specific Reaction: 1–3+ signal, trace background on the positive well and negative on the negative well.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Individual cells with diffuse cytoplasmic and membrane fluorescence and strong nuclear fluorescence.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with Bovine Adenovirus Type 1 (BAV-1) at 1–2+, but does not react with Bovine Adenovirus Type 5 (BAV-5), Bovine Coronavirus (BCV), Bovine Herpesvirus Type 1 (BHV-1/IBR), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Parainfluenza Virus Type 3 (PI-3), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

6 April 2023