

CERTIFICATE OF ANALYSIS

Bovine Adenovirus Type 5 (BAV-5)

FITC Conjugate


Catalog No.:	CJ-F-BAV5-10ML
Lot:	P190719-001
Expiration:	06 July 2031

Quality Control Method: Direct FA using BAV-5 2-well slide (catalog no. SLD-FAC-BAV5).

Specific Reaction: 3-4+ signal on the positive well and negative on the negative well, no background.

Other Comments: NA

QC Manager:



Jennifer Boren

Pattern Of Fluorescence: Diffuse and granular cytoplasmic fluorescence. Brightly fluorescing nuclear inclusions and nuclei. Some infected cells are multinucleated.

FA Cross-Reaction Screening: The raw material has also been screened by Direct FA and was found to react with both Bovine Adenovirus Type 1 (BAV-1) and Bovine Adenovirus Type 3 (BAV-3) at trace - 1+, but does not react with Bovine Coronavirus (BCV), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Infectious Bovine Rhinotracheitis (IBR/BHV-1), Bovine Parainfluenza virus type 3 (PI-3), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

12 July 2023