

CERTIFICATE OF ANALYSIS

Anti-Bovine IgG_{1,2} (H&L)


FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-BOVG-AP-10ML
Lot:	P211101-001
Expiration:	06 June 2026

Quality Control Method: Indirect FA using Bovine Respiratory Syncytial Virus (BRSV) 12-well slide (SLD-IFA-BRSV), positive control (PC-IFA-BRSV), and negative control (NC-IFA-BRSV).

Specific Reaction: 2–4+ signal, no background with the positive control at neat. Negative, trace background with the negative control.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: The pattern of fluorescence will vary depending on what bovine system is used.

FA Cross-Reaction Screening: The raw material has also been tested satisfactorily by VMRD's IFA systems for Bovine Adenovirus Type 1, 3, and 5 (BAV-1, 3, and 5), Bovine Coronavirus (BCV), Bovine Herpesvirus type 1 (BHV-1/IBR), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Caprine Arthritis Encephalitis Virus (CAEV), *Neospora caninum* (bovine origin), Parainfluenza Virus type 3 (PI-3), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

15 June 2022