

CERTIFICATE OF ANALYSIS

Bovine Viral Diarrhea Virus (BVDV)

FITC Conjugate

Catalog No.:	CJ-F-BVD-10ML
Lot:	P240301
Expiration:	02 April 2032

Quality Control Method: Direct FA using BVDV 2-well slide (catalog no. SLD-FAC-BVD).

Specific Reaction: 2–3+ signal, trace background on the positive well and negative, trace background on the negative well.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Bovine Adenovirus Type 1, 3, and 5 (BAV-1, 3, and 5), Bovine Coronavirus (BCV), Bovine Herpesvirus type 1 (BHV-1/IBR), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Parainfluenza Virus Type 3 (PI-3), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

12 April 2024