

CERTIFICATE OF ANALYSIS

## Anti-Canine IgG


FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-CANG-AP-10ML
Lot:	M230220-001
Expiration:	06 March 2031

**Quality Control Method:** Indirect FA using CDV 12-well slide (catalog no. SLD-IFA-CDV), and CDV IgG positive control (catalog no. PC-IFA-CDV-G).

**Specific Reaction:** 3–4+ signal with the positive control at neat, no background.

**Other Comments:** NA

**QC Manager:** 

David Hyatt

**Pattern Of Fluorescence:** The pattern of fluorescence will vary depending on what canine system is used.

**FA Cross-Reaction Screening:** The raw material has also been tested satisfactorily by VMRD's IFA systems for *Brucella canis*, Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), *Ehrlichia canis*, *Leishmania infantum*, *Borrelia burgdorferi* (Lyme Disease), *Neospora caninum* (canine origin), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Rickettsia rickettsii* (RMSF), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

23 May 2023