

CERTIFICATE OF ANALYSIS

Anti-Canine IgG

FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-CANG-AP-1ML
Lot:	P170725-001
Expiration:	17 August 2021

Quality Control Method: Indirect FA using CDV 12-well slides (catalog no. SLD-IFA-CDV), and CDV IgG positive control (catalog no. PC-IFA-CDV-G).

Specific Reaction: 2-4+ signal with the CDV positive control, no background. Negative with diluent control, no background.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 24 August 2017

Pattern of Fluorescence: The pattern of fluorescence will vary depending on what canine system was used.

FA Cross-Reaction Screening: This product and /or the raw material from which it was made has been tested in IFA and found to not react with Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus (CCV), *Brucella canis*, Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), *Ehrlichia canis*, *Borrelia burgdorferi* (Lyme Disease), *Leishmania infantum*, *Neospora caninum* (canine origin), *Rickettsia rickettsii* (RMSF), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR IN VITRO LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Anti-Immunoglobulin CJ\Canine IgG Affinity Purified\Canine IgG AP - 1 ml\CJ-F-CANG-AP-1ML_P170725-001_210817.docx

15 September 2017