

CERTIFICATE OF ANALYSIS

ANTI-CANINE IgM


FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-CANM-AP-10ML
Lot:	P210629-001
Expiration:	23 April 2025

Quality Control Method: Indirect FA using CDV 12-well slide (catalog no. SLD-IFA-CDV), and CDV IgM positive control (catalog no. PC-IFA-CDV-M).

Specific Reaction: 2-4+ signal with CDV IgM positive control at neat, trace background.
Negative with the diluent control, trace background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: The pattern of fluorescence will vary depending on what system is used.

FA Cross-Reaction Screening: The raw material has also been tested satisfactorily by VMRD's IFA systems for Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), *Ehrlichia canis*, *Leishmania infantum*, *Borrelia burgdorferi* (Lyme Disease), *Neospora caninum* (canine origin), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Rickettsia rickettsii* (RMSF), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

29 April 2021