

CERTIFICATE OF ANALYSIS

ANTI-CANINE IgM

FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-CANM-AP-1ML
Lot:	P161024-004
Expiration:	28 April 2021

Quality Control Method: Indirect FA using CDV 12-well slide (catalog no. SLD-IFA-CDV), and CDV IgM positive control (catalog no. PC-IFA-CDV-M).

Specific Reaction: 2-4+ signal with CDV IgM positive control at neat and no background.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 12 May 2017

Pattern of Fluorescence: The pattern of fluorescence will vary depending on what canine system was used.

FA Cross-Reaction Screening: The raw material has also been tested satisfactorily by VMRD's IFA systems for Anaplasma phagocytophilum, Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus (CCV), Brucella canis, Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), Ehrlichia canis, Borrelia burgdorferi (Lyme Disease), Leishmania infantum, Neospora caninum (canine origin), Rickettsia rickettsii (RMSF), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), Toxoplasma gondii, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Anti-Immunoglobulin CJ\Canine IgM\Canine IgM - 1 mL\CJ-F-CANM-AP-1ML_P161024-004_210428.docx 07 June 2021