

CERTIFICATE OF ANALYSIS

Canine Adenovirus (CAV)

FITC Conjugate

Catalog No.:	CJ-F-CAV-10ML
Lot:	P180809-001
Expiration:	22 July 2030

Quality Control Method: Direct FA using CAV-2 2-well slide (catalog no. SLD-FAC-CAV).

Specific Reaction: 1-3+ signal on the positive well and negative on the negative well. No

background.

Other Comments:

QC Manager:

David Hyatt

Pattern Of Fluorescence: Mostly large, single nuclear inclusion and some cytoplasmic fluorescence

in a few infected cells.

FA Cross-Reaction The raw material used in the production of this conjugate has also been Screening: screened by FA and has been found to react with Canine Adenovirus Type 1 (CAV-1) at a 1/200 dilution with a 2-4+ signal, trace background, but does not react with Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), Reovirus (REO), Rabies

Recombinant Nucleoprotein (rNP), Toxoplasma gondii, Vesicular Stomatitis

Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.