

## **CERTIFICATE OF ANALYSIS**

## Canine Coronavirus (CCV)

FITC Conjugate

Catalog No.:	CJ-F-CCV-10ML
Lot:	P210727-001
Expiration:	21 June 2025

Quality Control Method: Direct FA using CCV 2-well slide (catalog no SLD-FAC-CCV).

Specific Reaction: 2-4+ signal on the positive well and negative on the negative well, trace

background.

Other Comments: NA

QC Manager:

David Hyatt

Pattern Of Fluorescence: Large fluorescent syncytia with some single cells in between.

FA Cross-Reaction The raw material has also been screened by FA and does not react with

Screening: Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Distemper Virus

(CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus

Type 2 (CPI-2), Canine Parvovirus (CPV), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular

Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.