

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV)

FITC Conjugate

Catalog No.:	CJ-F-CDV-10ML
Lot:	P200505-001
Expiration:	18 March 2024

Quality Control Method: Direct FA on CDV 2-well slide (catalog no. SLD-FAC-CDV) and CDV negative blood smear (catalog no. SLD-BSN-CDV).

Specific Reaction: 3-4+ signal on the positive well and negative on the negative well. Negative on the negative blood smears.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane fluorescence.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with canine adenovirus type 1 and 2 (CAV 1 and 2), canine coronavirus (CCV), canine herpesvirus type 1 (CHV), canine parainfluenza virus type 2 (CPI-2), canine parvovirus (CPV), reovirus (REO), and Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, and vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

14 April 2020