

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV)

FITC Conjugate MAb

Catalog No.:	CJ-F-CDV-MAB-10ML
Lot:	P150727-001
Expiration:	6 August 2019

Quality Control Method: Direct FA on CDV 2-well slide (catalog no. SLD-FAC-CDV), CDV positive

blood smear (catalog no. SLD-BSP-CDV), and CDV negative blood smear

(catalog no. SLD-BSN-CDV).

Specific Reaction: 2-4+ fluorescence on the positive well and negative on the negative well,

no background. 2-4+ fluorescence on the positive smear and negative on

the negative smear, no background.

Other Comments: NA

QC Manager:

Will Thurman

QC Release Date: 7 August 2015

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane

fluorescence.

FA Cross-Reaction The raw material has also been screened by Indirect FA and does not react

Screening: with Canine Adenovirus Type 1 and 2 (CAV-1 and 2), Canine Coronavirus

(CCV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus

Type 2 (CPI-2), Canine Parvovirus (CPV), Reovirus (REO), Rabies

Recombinant Nucleoprotein (rNP), Toxoplasma gondii, Vesicular Stomatitis

Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1