

CERTIFICATE OF ANALYSIS

Canine Parainfluenza Virus Type 2 (CPI-2)

FITC Conjugate

Catalog No.:	CJ-F-CPI-10ML
Lot:	P170328-001
Expiration:	01 February 2025

Quality Control Method: Direct FA using CPI-2 2-well slide (catalog no. SLD-FAC-CPI).

Specific Reaction: 2-4+ signal on the positive well and negative on the negative well. No background.

Other Comments: NA

QC Manager:



David Hyatt

Pattern of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with canine herpesvirus type 1 (CHV-1) at trace to 1+ and *Toxoplasma gondii* at trace to 1+ with this ready to use conjugate, but does not react with canine adenovirus type 1 and 2 (CAV-1 and 2), canine coronavirus (CCV), canine distemper virus (CDV), canine parvovirus (CPV), reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Conjugate\CPI-2\CPI-2 - 10 ml\CJ-F-CPI-10ML_P170328-001_250201 (002).docx

04 February 2021