

CERTIFICATE OF ANALYSIS

Canine Parainfluenza Virus Type 2 (CPI-2)


FITC Conjugate

Catalog No.:	CJ-F-CPI-10ML
Lot:	P220525-001
Expiration:	22 June 2026

Quality Control Method: Direct FA using CPI-2 2-well slide (catalog no. SLD-FAC-CPI).

Specific Reaction: 3–4+ signal on the positive well and negative on the negative well. No background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with *Toxoplasma gondii* at Trace-1+, but does not react with canine adenovirus type 1 and 2 (CAV-1 and 2), canine coronavirus (CCV), canine distemper virus (CDV), canine herpesvirus type 1 (CHV-1), canine parvovirus (CPV), Rabies Recombinant Nucleoprotein (rNP), reovirus (REO), vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

11 July 2022