

## **CERTIFICATE OF ANALYSIS**

## Clostridium sordellii

FITC Conjugate

Catalog No.:	CJ-F-CSO-10ML
Lot:	P200121-001
Expiration:	03 January 2024

Quality Control Method: Direct FA using Clostridium sordellii 12-well slide (catalog no. SLD-IFA-

CSO).

Specific Reaction: 1-3+ signal with no background.

Other Comments: NA

QC Manager:

David Hyatt

Pattern Of Fluorescence: Individual rod-shaped bacteria. Some clumped bacteria. Flagella

fluorescence alone is not indicative of a positive signal.

FA Cross-Reaction The raw material has also been screened by FA and does not react with

Screening: C. novyi, C. septicum, or C. chauvoei. There was trace reaction with

C. perfringens.

Technical Data Sheet Version: 2