

CERTIFICATE OF ANALYSIS

Anti-Equine IgG

FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-EQUG-AP-10ML
Lot:	P201215-001
Expiration:	12 February 2025

Quality Control Method: Indirect FA on Equine Herpesvirus Type 1 (EHV-1/ ERV) 12-Well Slide (catalog no. SLD-IFA-ERV), and positive control (catalog no. PC-IFA-ERV-EQ).

Specific Reaction: 2–4+ signal with the positive control at neat, trace–1+ background. Negative with the diluent control, no background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: The pattern of fluorescence will vary depending on what equine system was used.

FA Cross-Reaction Screening: The raw material has also been tested satisfactorily by VMRD's IFA systems for *Anaplasma phagocytophilum*, Equine Arteritis Virus (EAV), Equine Herpesvirus type 1 (EHV-1/ERV), *Borrelia burgdorferi* (Lyme disease), *Rickettsia rickettsii* (RMSF), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV)

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

24 February 2021