

CERTIFICATE OF ANALYSIS

Equine Herpesvirus Type 1 (EHV-1/ERV)

FITC Conjugate

Catalog No.:	CJ-F-ERV-10ML
Lot:	P170316-001
Expiration:	11 April 2021

Quality Control Method: Direct FA using EHV-1/ERV 2-well slide (catalog no. SLD-FAC-ERV).

Specific Reaction: 3-4+ signal on the positive well, negative on the negative well with no background.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 11 April 2017

Pattern of Fluorescence: Strong undifferentiated cytoplasmic fluorescence. Single cells with occasional nuclear inclusions and some small syncytia. Some intense fluorescence of dying rounded cells.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Equine Arteritis Virus (EAV), Equine Infectious Anemia Virus (EIAV), reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Conjugate\ERV\ERV - 10 ml\CJ-F-ERV-10ML_P170316001_210411.docx

18 April 2017