

CERTIFICATE OF ANALYSIS

## Anti-Feline IgM (heavy chain specific)

FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-FELM-AP-10ML
Lot:	P160115-003
Expiration:	10 February 2020

**Quality Control Method:** Indirect FA using *Bartonella henselae* 12-well slide (catalog no. SLD-IFA-BH) and IgM positive control (catalog no. PC-IFA-BH-M).

**Specific Reaction:** 3-4+ signal with the positive control at neat, no background. Negative with the diluent control.

**Other Comments:** NA

**QC Manager:**  Will Thurman

**QC Release Date:** 11 February 2016

**Pattern Of Fluorescence:** The pattern of fluorescence will vary depending on what feline system was used.

**FA Cross-Reaction Screening:** The raw material has also been tested satisfactorily by VMRD's IFA systems for *Bartonella Henselae*, Feline Calicivirus (FCV), Feline Leukemia Virus (FeLV), Feline Herpesvirus (FHV-1/FVR), Feline Immunodeficiency Virus (FIV), Feline Infectious Peritonitis type 1 and 2 (FIP-1 and 2), Feline Panleukopenia Virus (FPLV), Reovirus (REO), and Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version #1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Anti-Immunoglobulin CJ\Feline IgM\Feline IgM - 10 ml\CJ-F-FELM-AP-10ML\_P160115-003\_200210.docx

15 February 2016