

CERTIFICATE OF ANALYSIS

Anti-Feline IgM (heavy chain specific)


FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-FELM-AP-10ML
Lot:	P220530-001
Expiration:	17 August 2030

Quality Control Method: Indirect FA using *Bartonella henselae* 12-well slide (catalog no. SLD-IFA-BH) and IgM positive control (catalog no. PC-IFA-BH-M).

Specific Reaction: 4+ signal with the positive control at neat. Negative with the diluent control. No background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: The pattern of fluorescence will vary depending on what system is used.

FA Cross-Reaction Screening: The raw material has also been tested satisfactorily by VMRD's IFA systems for *Bartonella Henselae*, Feline Calicivirus (FCV), Feline Herpesvirus (FHV-1/FVR), Feline Infectious Peritonitis type 1 and 2 (FIP-1 and 2), Feline Panleukopenia Virus (FPLV), Feline Immunodeficiency Virus (FIV), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

24 August 2022