

CERTIFICATE OF ANALYSIS

Feline Infectious Peritonitis Virus (FIP)

FITC Conjugate

Catalog No.:	CJ-F-FIP-10ML
Lot:	P210420-001
Expiration:	16 June 2025

Quality Control Method: Direct FA using FIP-1 2-well slide (catalog no. SLD-FAC-FIP1) and FIP-2 2-well slide (catalog no. SLD-FAC-FIP2).

Specific Reaction: FIP-1 had 2–4+ signal on the positive well and was negative on the negative well, both with trace background. FIP-2 had trace–3+ signal on the positive well and was negative on the negative well with no background.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Multiple syncytia, plaques and some individual cells with granular cytoplasmic fluorescence.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with *Toxoplasma gondii* with 1+ signal, but does not react with feline calicivirus (FCV), feline leukemia virus (FeLV), feline herpesvirus (FHV-1/FVR), feline panleukopenia virus (FPLV), feline immunodeficiency virus (FIV), rabies recombinant nucleoprotein (rNP), reovirus (REO), vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

8 July 2021