

CERTIFICATE OF ANALYSIS

Anti-Murine IgG


FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-MURG-AP-10ML
Lot:	P200122-001
Expiration:	29 September 2031

Quality Control Method: Indirect FA using CDV 12-well slide (catalog no. SLD-IFA-CDV), CDV MAb CDV-NP (catalog no. CDV-NP), and isotype control IgG_{2b}.

Specific Reaction: 2-3+ signal with the CDV-NP monoclonal at 10 µg/ml, no background. Negative with the isotype control control, no background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: The pattern of fluorescence will vary depending on what system was used.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with *Toxoplasma gondii* at trace-1+.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

18 October 2023