

CERTIFICATE OF ANALYSIS

Porcine Circovirus Type 1 & 2 (PCV1&2)

FITC Conjugate

Catalog No.:	CJ-F-PCV1&2-10ML
Lot:	P200902-001
Expiration:	23 September 2024

Quality Control Method: Direct FA using PCV-1 2-well slide (catalog no. SLD-FAC-PCV1) and PCV-2

2-well slide (catalog no. SLD-FAC-PCV2).

Specific Reaction: 2–4+ fluorescence on the positive wells of the PCV-1 and PCV-2 slides.

Negative with trace background on the negative wells of the PCV-1 slides, and negative with no background on the negative wells of the PCV-2

slides.

Other Comments:

NA

QC Manager:

David Hyatt

Pattern Of Fluorescence: Intense whole nuclear and variable sized but discrete cytoplasmic

inclusions.

FA Cross-Reaction The raw material has also been screened by FA and does not react with

Screening: porcine adenovirus (PAV), porcine hemagglutinating encephalomyelitis

virus (PHEV), porcine parvovirus (PPV), porcine reproductive and

respiratory syndrome virus (PRRSV) North American and European strains,

transmissible gastroenteritis virus (TGEV), rabies recombinant nucleoprotein (rNP), reovirus (REO), *Toxoplasma gondii*, vesicular

stomatitis virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: 4

FOR IN VITRO LABORATORY USE ONLY.

1 October 2020