

CERTIFICATE OF ANALYSIS

Porcine Circovirus Type 2 (PCV-2)

FITC Conjugate

Catalog No.:	CJ-F-PCV2-10ML
Lot:	P151026-001
Expiration:	27 October 2019

Quality Control Method: Direct FA using PCV-2 2-well slide (catalog no. SLD-FAC-PCV2).

Specific Reaction: 2-4+ signal on the positive well, negative on the negative well with no

background noted on either well.

Other Comments: NA

Will Thurman QC Manager:

QC Release Date: 27 October 2015

Pattern Of Fluorescence: Intense whole nuclear and variable sized but discrete cytoplasmic

inclusions.

FA Cross-Reaction The raw material has also been screened by FA and does not react with

Porcine Adenovirus (PAV), Porcine Cytomegalovirus (PCMV), Porcine Screening:

Circovirus Type 1 (PCV-1), Porcine Hemagglutinating Encephalomyelitis

Virus (PHEV), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome Virus (PRRSV) North American and European

strains, Transmissible Gastroenteritis Virus (TGEV), Reovirus (REO), Rabies

Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus (VSV) Indiana

and New Jersey strains.

Technical Data Sheet Version: Version # 1