

CERTIFICATE OF ANALYSIS

Anti-Porcine IgG (H&L)

FITC Conjugate Affinity Purified

Catalog No.:	CJ-F-PORG-AP-10ML
Lot:	P160115-004
Expiration:	29 January 2020

Quality Control Method: Indirect FA using PCV-2 12-well slide (catalog no. SLD-IFA-PCV2), positive control (catalog no. PC-IFA-PCV2), and negative control (catalog no. NC-IFA-PCV2).

Specific Reaction: 2-4+ signal with PCV2 positive control at neat, no background. 1+ endpoint reached in three two-fold dilutions.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 01 February 2016

Pattern Of Fluorescence: The pattern of fluorescence will vary depending on what porcine system was used.

FA Cross-Reaction Screening: The raw material has also been tested satisfactorily by VMRD's IFA systems for Porcine Adenovirus (PAV), Porcine Circovirus Type 2 (PCV-2), Porcine Hemagglutinating Encephalomyelitis Virus (PHEV), Porcine Parvovirus (PPV), Porcine Reproductive and Respiratory Syndrome virus (PRRSV) North American and European strains, Vesicular Stomatitis Virus (VSV) Indiana and New Jersey strains and Transmissible Gastroenteritis Virus (TGEV). This FITC conjugate reacts weakly (trace to 1+) with *Toxoplasma gondii*.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Anti-Immunoglobulin CJ\Pig IgG\Pig IgG - 10 ml\CJ-F-PORG-AP-10ML_P160115-004_200129.docx

03 February 2016