

CERTIFICATE OF ANALYSIS

Porcine Parvovirus (PPV)


FITC Conjugate

Catalog No.:	CJ-F-PPV-10ML
Lot:	P231206
Expiration:	12 December 2031

Quality Control Method: Direct FA using PPV 2-well slide (catalog no. SLD-FAC-PPV).

Specific Reaction: 2-4+ signal on the positive well and negative on the negative well. No background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Fluorescence limited almost entirely to nucleus of cell. Some granular cytoplasmic staining.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with porcine cytomegalovirus (PCMV) 3+ at 1/50, but does not react with porcine adenovirus (PAV), porcine circovirus type 1 and 2 (PCV-1 and 2), porcine hemagglutinating encephalomyelitis virus (PHEV), porcine reproductive and respiratory syndrome virus North American and European strains (PRRSV), transmissible gastroenteritis virus (TGEV), rabies Recombinant Nucleoprotein (rNP), reovirus (REO), *Toxoplasma gondii*, vesicular stomatitis virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

18 December 2023