

**CERTIFICATE OF ANALYSIS**

## Reovirus (REO)

FITC Conjugate

Catalog No.:	CJ-F-REO-10ML
Lot:	P180219-002
Expiration:	30 June 2026

**Quality Control Method:** Direct FA using REO 2-well slide (catalog no. SLD-FAC-REO).

**Specific Reaction:** 3–4+ signal on the positive well and negative on the negative well. No background.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Individual cells or plaques with diffuse cytoplasmic fluorescence and characteristic bright cytoplasmic inclusions ranging in size from medium to large.

**FA Cross-Reaction Screening:** The raw material has also been screened by FA and does not react with Bovine Adenovirus Type 1, 3, and 5 (BAV-1, 3, and 5), Bovine Coronavirus (BCV), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Bovine Herpesvirus type 1 (BHV-1/IBR), Parainfluenza Virus Type 3 (PI-3), Rabies Recombinant Nucleoprotein (rNP), Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Master Documents\COA (VMRD)\Conjugate\Reovirus\Reovirus - 10 ml\CJ-F-REO-10ML\_P180219-002\_260630.docx

01 July 2022