

CERTIFICATE OF ANALYSIS

## Reovirus (REO)

FITC Conjugate

Catalog No.:	CJ-F-REO-10ML
Lot:	P220209-001
Expiration:	09 February 2026

Quality Control Method: Direct FA using REO 2-well slide (catalog no. SLD-FAC-REO).

Specific Reaction: 3–4+ signal on the positive well and negative with trace background on the negative well.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells or plaques with diffuse cytoplasmic fluorescence and characteristic bright cytoplasmic inclusions ranging in size from medium to large.

FA Cross-Reaction Screening: The raw material has also been screened by FA and does not react with Bovine Adenovirus Type 1, 3, and 5 (BAV-1, 3, and 5), Bovine Coronavirus (BCV), Bovine Herpesvirus type 1 (BHV-1/IBR), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Parainfluenza Virus Type 3 (PI-3), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

25 February 2022