

CERTIFICATE OF ANALYSIS

Clostridium novyi

FITC Conjugate

Catalog No.:	CJ-F-CNO-10ML
Lot:	P180502-005
Expiration:	15 June 2026

Quality Control Method: Direct FA using the C. novyi well of a Clostridium spp. 4-way slide. (catalog

no. SLD-IFA-C4).

Specific Reaction: 2-3+ signal with no background.

Other Comments: NA

QC Manager:

David Hyatt

David Wyath

Pattern Of Fluorescence: Individual rod-shaped bacteria. Some clumped bacteria. Flagella

fluorescence alone is not indicative of a positive signal.

FA Cross-Reaction The raw material has also been screened by FA and does not react with C.

Screening: chauvoei, C. perfringens, C. septicum, or C. sordellii.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.