

CERTIFICATE OF ANALYSIS

FA Conjugate Diluting Buffer

Catalog No.:	FACDB-100ML
Lot:	P160829-005
Expiration:	31 October 2020

Quality Control Method: Direct FA using 1/50 dilution of 50X Conjugates and the corresponding slides on each of the following VMRD systems: Bovine PI-3 (catalog no. CJ-F-PI3-50X, and SLD-IFA-PI3), Canine CDV (catalog no. CJ-F-CDV-50X, and SLD-IFA-CDV), Equine/Llama ERV (catalog no. CJ-F-ERV-50X, and SLD-IFA-ERV), Feline FVR (catalog no. CJ-F-FVR-50X, and SLD-IFA-FVR), and Porcine TGE (catalog no. CJ-F-TGE-50X, and SLD-IFA-TGE).

Specific Reaction: 2-4+ signal for all systems at 1/50 with no additional background compared to the previous lot. CDB alone produced no background and no signal with all systems.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 01 November 2016

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Buffers\Conjugate Diluting Buffer\100 ml\FACDB-100ML_P160829-005_201031.docx

03 November 2016