

CERTIFICATE OF ANALYSIS


FA Conjugate Diluting Buffer

Catalog No.:	FACDB-100ML
Lot:	P201014-001
Expiration:	07 October 2024

Quality Control Method: Direct FA using 1/50 dilution of 50X conjugates and the corresponding slides on each of the following VMRD systems: Bovine BRSV (catalog no. CJ-F-BRSV-50X, and BRSV in-house slide), Canine CDV (catalog no. CJ-F-CDV-50X, and SLD-IFA-CDV), Equine/Llama ERV (catalog no. CJ-F-ERV-50X, and SLD-IFA-ERV), Feline FVR (catalog no. CJ-F-FVR-50X, and SLD-IFA-FVR), and Porcine TGE (catalog no. CJ-F-TGE-50X, and SLD-IFA-TGE).

Specific Reaction: Each system produced in the range of 2–4+ fluorescence with 0–trace background at 1/50. The diluent control was negative with no background on all systems except porcine TGE, which produced trace fluorescence.

Other Comments: NA

QC Manager: 

David Hyatt

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

15 October 2020