

CERTIFICATE OF ANALYSIS

## 4X FA Powdered Rinse Buffer

Catalog No.:	FARB-4X
Lot:	P180921-001
Expiration:	05 October 2022

**Quality Control Method:** Reconstitute in 1 liter of deionized water to make 4X concentrate. Using pH meter, determine the pH of the solution. Solution should be between pH 9.0 and pH 9.5. Dilute to 1X. Direct FA using FIP conjugate (catalog no. CJ-F-FIP-1ML or 10ML) and FIP2 2-well control slide (catalog no. SLD-FAC-FIP2).

**Specific Reaction:** Reconstituted 4X and 1X rinse buffer was clear. QC rinse buffer produced results similar to the reference rinse buffer.

pH: 9.30

**Other Comments:** NA

**QC Manager:** Chris Cay



**QC Release Date:** 08 October 2018

**Technical Data Sheet Version:** Version # 1

FOR IN VITRO LABORATORY USE ONLY.

17 October 2018