

CERTIFICATE OF ANALYSIS

4X FA Powdered Rinse Buffer

Catalog No.:	FARB-4X
Lot:	P201023-001
Expiration:	16 November 2024

Quality Control Method: Reconstitute in 1 liter of deionized water to make 4X concentrate. Using pH meter, determine the pH of the solution. Solution should be between pH 9.0 and pH 9.5. Dilute to 1X. Direct FA using FIP conjugate (catalog no. CJ-F-FIP-10ML) and FIP2 2-well control slide (catalog no. SLD-FAC-FIP2).

Specific Reaction: Reconstituted 4X and 1X rinse buffer was clear. QC rinse buffer produced results similar to reference rinse buffer.

pH: 9.10

Other Comments: NA

QC Manager:



David Hyatt

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

23 November 2020