

CERTIFICATE OF ANALYSIS

4X FA Powdered Rinse Buffer

Catalog No.:	FARB-4X
Lot:	P240212
Expiration:	19 March 2032

Quality Control Method: Reconstituted in 1 liter of deionized water to make 4X concentrate and confirmed pH of reconstituted solution to be between 9.0 and 9.5. Diluted 4X solution to 1X and used for indirect FA using FIP-1 slides (catalog no. SLD-IFA-FIP1), FIP-1 positive control (catalog no. PC-IFA-FIP1), negative control (catalog no. NC-IFA-FIP1), and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: QC rinse buffer produced results similar to reference rinse buffer.

pH: 9.02

Other Comments: NA

QC Manager:



David Hyatt

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

28 March 2024