

**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vtech@vmrd.comWeb site: <http://www.vmrd.com>

Certificate of Analysis

Bluetongue Virus (BTV)

Negative Control Serum

CATALOG NO.: NC-IFA-BTV**VOLUME:** 1 ml**LOT:** V98010**EXPIRATION:** 23 October 2024**VIRUS:** Bluetongue Virus (BTV)**DESCRIPTION:** Bovine serum diluted in PBS, 1% BSA, 0.09% sodium azide.**INTENDED USE:** As negative control serum in detection of antibody to Bluetongue Virus (BTV) by indirect FA technique. This serum should be used undiluted to demonstrate nonspecific fluorescence, if any.**QUALITY CONTROL METHOD:** Indirect FA using BTV 12-well slide (catalog no. SLD-IFA-BTV), positive control (catalog no. PC-IFA-BTV), and anti-bovine IgG1,2 affinity purified FITC conjugate (catalog no. CJ-F-BOVG-AP-10).**Specific Reaction:** Negative using the negative control with trace-1+ background.**Other Reactions or Comments:** NA**PATTERN OF FLUORESCENCE:** None with our system.**STORAGE:** Store at 2-7°C. DO NOT FREEZE!

FOR *IN VITRO* LABORATORY USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.