

CERTIFICATE OF ANALYSIS

Brucella canis

FA Negative Control

Catalog No.:	NC-IFA-CB
Lot:	P240407
Expiration:	18 April 2032

Quality Control Method: Indirect FA using *Brucella canis* 12-well slides (catalog no. SLD-IFA-CB), positive control (catalog no. PC-IFA-CB), and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: Negative with no background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: None with our system.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

22 April 2024