

CERTIFICATE OF ANALYSIS

*Neospora caninum*

FA Negative Control

Catalog No.:	NC-IFA-NC-CAN
Lot:	P181220-001
Expiration:	01 April 2023

**Quality Control Method:** Indirect FA using *Neospora caninum* 12-well slide (catalog no. SLD-IFA-NC), canine positive control (catalog no. PC-IFA-NC-CAN), and Anti-Canine FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

**Specific Reaction:** Negative with trace to 1+ apical staining background.

**Other Comments:** When using *Neospora caninum* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

**QC Manager:** Chris Cay



**QC Release Date:** 01 April 2019

**Pattern Of Fluorescence:** None with our system.

**Technical Data Sheet Version:** Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

22 April 2019