

CERTIFICATE OF ANALYSIS

Toxoplasma gondii

FA Negative Control

Catalog No.:	NC-IFA-TOXO-FEL
Lot:	P170313-001
Expiration:	31 March 2021

Quality Control Method: Indirect FA using *Toxoplasma gondii* 12-well slide (catalog no. SLD-IFA-TOXO), IgG positive control (catalog no. PC-IFA-TOXO-FEL-G) and anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-1ML or 10ML).

Specific Reaction: Negative with trace-1+ apical signal.

Other Comments: When using *Toxoplasma gondii* systems it is recommended to use FASDB as a diluent.

QC Manager:  Will Thurman

QC Release Date: 31 March 2017

Pattern of Fluorescence: None with our system.

Technical Data Sheet Version: Version #2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\Negative Controls\Toxoplasma\Feline\NC-IFA-TOXO-FEL_P170313-001_200331.docx

24 April 2017