

CERTIFICATE OF ANALYSIS

*Toxoplasma gondii*

FA Negative Control

|              |                 |
|--------------|-----------------|
| Catalog No.: | NC-IFA-TOXO-FEL |
| Lot:         | P191104-002     |
| Expiration:  | 04 April 2031   |

**Quality Control Method:** Indirect FA using *Toxoplasma gondii* 12-well slide (catalog no. SLD-IFA-TOXO), IgG positive control (catalog no. PC-IFA-TOXO-FEL-G) and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** Negative with 1+ apical signal.

**Other Comments:** When using *Toxoplasma gondii* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

**QC Manager:** 

David Hyatt

**Pattern Of Fluorescence:** None with our system.

Technical Data Sheet Version: 2