

CERTIFICATE OF ANALYSIS

Toxoplasma gondii

FA Negative Control

Catalog No.:	NC-IFA-TOXO-FEL
Lot:	P231002
Expiration:	06 October 2031

Quality Control Method: Indirect FA using *Toxoplasma gondii* 12-well slide (catalog no. SLD-IFA-TOXO), IgG positive control (catalog no. PC-IFA-TOXO-FEL-G) and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: Negative with trace-1+ apical signal.

Other Comments: When using *Toxoplasma gondii* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

QC Manager: 

David Hyatt

Pattern Of Fluorescence: None with our system.

Technical Data Sheet Version: 2