

CERTIFICATE OF ANALYSIS

Porcine Circovirus Type 2 (PCV-2)


Antiserum

Catalog No.:	PAB-PCV2
Lot:	P150204-001
Expiration:	12 December 2030

Quality Control Method: Indirect FA using PCV-2 12-well slide (catalog no. SLD-IFA-PCV2) and anti-porcine IgG affinity purified FITC conjugate (catalog no. CJ-F-PORG-AP-10ML).

Specific Reaction: 2-4+ signal, trace background at 1/50 with an endpoint greater than 1/400.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Granular and diffuse nuclear fluorescence with some nuclear inclusions and occasional smooth and undifferentiated cytoplasmic fluorescence.

FA Cross-Reaction Screening: The Antiserum has also been screened by indirect FA and has been found to react with porcine adenovirus (PAV) trace to 1+ at 1/100, but does not react with porcine circovirus type 1 (PCV-1), porcine hemagglutinating encephalomyelitis virus (PHEV), porcine parvovirus (PPV), porcine reproductive and respiratory syndrome virus (PRRSV) North American and European strains, transmissible gastroenteritis virus (TGEV), reovirus (REO), rabies Recombinant Nucleoprotein (rNP), vesicular stomatitis virus (VSV) Indiana and New Jersey strains.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Antisera\PCV\PAB-PCV2_P150204-001_301212.docx

29 December 2022