

CERTIFICATE OF ANALYSIS

Parainfluenza Virus Type 3 (PI-3)


Antiserum

Catalog No.:	PAB-PI3
Lot:	P171017-003
Expiration:	07 October 2025

Quality Control Method: Indirect FA using PI-3 2-well slide (catalog no. SLD-FAC-PI3), and anti-caprine IgG FITC conjugate (catalog no. CJ-F-CAPG-10ML).

Specific Reaction: 2–4+ signal at a 1/400 dilution with no background.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Small, bright spherical cytoplasmic inclusions in syncytia and single cells.

FA Cross-Reaction Screening: The Antiserum has also been screened by Indirect FA and does not react with Bovine Adenovirus Type 1, 3, and 5 (BAV-1, 3, and 5), Bovine Coronavirus (BCV), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), Bovine Herpesvirus type 1 (BHV-1/IBR), Reovirus (REO), Rabies Recombinant Nucleoprotein (rNP), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: Version 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Antisera\PI-3\PAB-PI3_P171017-003_251007.docx

03 November 2021