

CERTIFICATE OF ANALYSIS

# Parainfluenza Virus Type 3 (PI-3)


Antiserum

Catalog No.:	PAB-PI3
Lot:	P211009-001
Expiration:	04 November 2025

**Quality Control Method:** Indirect FA using PI-3 2-well slide (catalog no. SLD-FAC-PI3), and anti-caprine IgG FITC conjugate (catalog no. CJ-F-CAPG-10ML).

**Specific Reaction:** 1–3+ signal with the antisera at 1/400 dilution, no background.

**Other Comments:** NA

**QC Manager:** 

David Hyatt

**Pattern Of Fluorescence:** Small, bright spherical cytoplasmic inclusions in syncytia and single cells.

**FA Cross-Reaction Screening:** The Antiserum has also been screened by FA and does not react with Bovine Adenovirus Type 1, 3, and 5 (BAV-1, 3, and 5), Bovine Coronavirus (BCV), Bovine Herpesvirus type 1 (BHV-1/IBR), Bovine Leukemia Virus (BLV), Bovine Parvovirus (BPV), Bovine Respiratory Syncytial Virus (BRSV), Bluetongue Virus (BTV), Bovine Viral Diarrhea Virus (BVDV), *Neospora caninum*, Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), *Toxoplasma gondii*, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

24 October 2023